

Product Datasheet

Mouse Monoclonal anti-human CD32 b and c (orb1153235)

Description	Mouse Monoclonal anti-human CD32 b and c
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human CD32 b/c
Storage	Reconstitute the antibody with 500 ul sterile PBS and the final concentration is 200ug/ml. Reconstituted antibody can also be aliquoted and stored frozen at -20 for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG2
Clonality	Monoclonal
Dilution Range	ELISA (1 to 1000-2000)
Expiration Date	12 months from date of receipt.