

Product Datasheet

Mouse Monoclonal anti-human Dkk-4 (orb1153242)

Description	Mouse Monoclonal anti-human Dkk-4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human Dkk-4
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 500 ug/ml. Reconstituted antibody can also be aliquoted and stored frozen at -20 for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG2
Clonality	Monoclonal
Dilution Range	WB (1 to 500-1000)
Expiration Date	12 months from date of receipt.