

Product Datasheet

Mouse Monoclonal anti-human KIR2DL4 (orb1153272)

Description	Mouse Monoclonal anti-human KIR2DL4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human recombinant KIR2DL4
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 500ug/ml. Reconstituted antibody can also be aliquoted and stored frozen at -20 for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles
Note	For research use only
Isotype	Mouse IgG2
Clonality	Monoclonal
Dilution Range	FC (1 to 50-100)
Expiration Date	12 months from date of receipt.