

Product Datasheet

Mouse Monoclonal anti-human IL-3R alpha EC (CD123) (orb1153312)

Description	Mouse Monoclonal anti-human IL-3R alpha EC (CD123)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant IL-3RalphaEC domain
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 500ug/ml. Reconstituted antibody can also be aliquoted and stored frozen at -20 for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Dilution Range	WB (1 to 500-1000), IHC-Fr (1 to 20-200)
Expiration Date	12 months from date of receipt.