

## Product Datasheet

### Mouse Monoclonal anti-human FGF-19 (orb1153453)

|                            |                                                                                                                                                                                                                                                                      |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse Monoclonal anti-human FGF-19                                                                                                                                                                                                                                   |
| <b>Species/Host</b>        | Mouse                                                                                                                                                                                                                                                                |
| <b>Reactivity</b>          | Human                                                                                                                                                                                                                                                                |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                                                         |
| <b>Tested Applications</b> | IHC, WB                                                                                                                                                                                                                                                              |
| <b>Immunogen</b>           | Human recombinant FGF-19                                                                                                                                                                                                                                             |
| <b>Storage</b>             | Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 500ug/ml. Reconstituted antibody can also be aliquoted and stored frozen at -20 for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                |
| <b>Isotype</b>             | Mouse IgG1                                                                                                                                                                                                                                                           |
| <b>Clonality</b>           | Monoclonal                                                                                                                                                                                                                                                           |
| <b>Dilution Range</b>      | WB (1 to 500-1000), IHC-P (1 to 50-100)                                                                                                                                                                                                                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                                      |