

Product Datasheet

Mouse Monoclonal anti-human FGF-17 (orb1153454)

Description	Mouse Monoclonal anti-human FGF-17
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Human recombinant FGF-17
Storage	Reconstitute the antibody with 200ul sterile PBS and the final concentration is 500ug/ml. Reconstituted antibody can also be aliquoted and stored frozen at -20 for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG2
Clonality	Monoclonal
Dilution Range	WB (1 to 500-1000), IHC-Fr (1 to 50-200)
Expiration Date	12 months from date of receipt.