

Product Datasheet

Mouse Monoclonal anti-human FGF-10 (orb1153457)

Description	Mouse Monoclonal anti-human FGF-10
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Immunogen	Human recombinant FGF-10
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 500ug/ml. Reconstituted antibody can also be aliquoted and stored frozen at -20 for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG3
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.