

Product Datasheet

Mouse Monoclonal anti-Human CD80 (B7-1) (orb1153689)

Description	Mouse Monoclonal anti-Human CD80 (B7-1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Human recombinant CD80 EC domain
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 500 ug/ml. Reconstituted antibody can also be aliquoted and stored frozen at -20 for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Dilution Range	WB (1 to 1000), IHC (1 to 50-100)
Expiration Date	12 months from date of receipt.