

Product Datasheet

Mouse Monoclonal anti-Human Aggrecan (G1-G2 domain) (orb1153707)

Description	Mouse Monoclonal anti-Human Aggrecan (G1-G2 domain)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human Aggrecan (G1-G2 domain)
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 500 ug/ml. The reconstituted antibody can also be aliquoted and stored frozen at -20 for at least six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG2
Clonality	Monoclonal
Dilution Range	WB (1 to 500-1000), IHC-P (1 to 40-100)
Expiration Date	12 months from date of receipt.