

Product Datasheet

Mouse Monoclonal anti-Human VCAM-1 (CD54) (orb1153739)

Description	Mouse Monoclonal anti-Human VCAM-1 (CD54)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Immunogen	LPS activated HUVECs
Storage	Reconstitute the antibody with 500 ul sterile PBS and the final concentration is 200ug/ml. Reconstituted antibody can also be aliquoted and stored frozen at -20 for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.