

Product Datasheet

Mouse Monoclonal anti-Rat CD11a Antibody (orb1153855)

Description	Mouse Monoclonal anti-Rat CD11a Antibody
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP
Immunogen	PHA-stimulated rat splenocytes and rat thymic lymphoma FTL-43
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.