

Product Datasheet

Mouse Monoclonal anti-RANTES (CCL5) Antibody (orb1153857)

Description	Mouse Monoclonal anti-RANTES (CCL5) Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant RANTES (CCL5)
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.