

Product Datasheet

Mouse Monoclonal anti-Pan-Shank Antibody (orb1153859)

Description	Mouse Monoclonal anti-Pan-Shank Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Fusion protein corresponding to amino acids 84-309 (SH3/PDZ domains) of rat Shank2
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.