

Product Datasheet

Rat Monoclonal anti-NGAL (Lipocalin-2) Antibody (orb1153860)

Description	Rat Monoclonal anti-NGAL (Lipocalin-2) Antibody
Species/Host	Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant human NGAL
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rat IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.