

Product Datasheet

Rat Monoclonal anti-Mouse/Human PNAd Antibody (orb1153873)

Description	Rat Monoclonal anti-Mouse/Human PNAd Antibody
Species/Host	Rat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP, WB
Immunogen	Mouse lymph node stromal cells
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rat IgM
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.