

Product Datasheet

Mouse Monoclonal anti-Mouse/Human NGF Antibody (orb1153874)

Description	Mouse Monoclonal anti-Mouse/Human NGF Antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	recombinant human NGF (a.a. 46-62)
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.