

Product Datasheet

Rat Monoclonal anti-Mouse/Human Ki-67 Antibody (orb1153876)

Description	Rat Monoclonal anti-Mouse/Human Ki-67 Antibody
Species/Host	Rat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	E. coli expressed, N-terminal His-Thioredoxin-tagged, partial mKi-67 (1816-2163 aa) recombinant protein.
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rat IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.