

Product Datasheet

Rat Monoclonal anti-Mouse/Human IL-5 Antibody (orb1153877)

Description	Rat Monoclonal anti-Mouse/Human IL-5 Antibody
Species/Host	Rat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, WB
Immunogen	Partially-purified T cell clone supernatant
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rat IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.