

Product Datasheet

Rat Monoclonal anti-Mouse/Human GL7 Antigen (T and B cell Activation Marker) Antibody (orb1153878)

Description	Rat Monoclonal anti-Mouse/Human GL7 Antigen (T and B cell Activation Marker) Antibody
Species/Host	Rat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, IP
Immunogen	LPS activated DBA/J mouse B cells
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rat IgM
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.