

Product Datasheet

Mouse Monoclonal anti-Mouse/Human Galectin-3 (Mac-2) Antibody (orb1153879)

Description	Mouse Monoclonal anti-Mouse/Human Galectin-3 (Mac-2) Antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human Recombinant (partial), amino acids 151-251
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.