

Product Datasheet

Rat Recombinant anti-Mouse/Human CD44 Recombinant Antibody (orb1153881)

Description	Rat Recombinant anti-Mouse/Human CD44 Recombinant Antibody
Species/Host	Rat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Recombinant mouse/human CD44
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rat IgG2a
Clonality	Recombinant
Expiration Date	12 months from date of receipt.