

Product Datasheet

Mouse Monoclonal anti-Mouse TCR V α 4 Antibody (orb1153904)

Description	Mouse Monoclonal anti-Mouse TCR V α 4 Antibody
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	C57BL/6 thymocytes expanded in vitro in anti-TCR α d antibody coated plates
Storage	Reconstitute the antibody with 200 μ l sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.