

Product Datasheet

Armenian Hamster Monoclonal anti-Mouse TCR v β 8.3 Antibody (orb1153906)

Description	Armenian Hamster Monoclonal anti-Mouse TCR v β 8.3 Antibody
Species/Host	Hamster
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	3.L2 T Cell Clone Expressing the v β 8.3 TCR
Storage	Reconstitute the antibody with 200 μ l sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Armenian Hamster IgG
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.