

Product Datasheet

Rat Monoclonal anti-Mouse TCR V β 8.1, 8.2 Antibody (orb1153907)

Description	Rat Monoclonal anti-Mouse TCR V β 8.1, 8.2 Antibody
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	TCR receptor from T cell hybridoma DO-11.10
Storage	Reconstitute the antibody with 200 μ l sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rat IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.