

Product Datasheet

Mouse Monoclonal anti-Mouse Qa-2 Antibody (orb1153909)

Description	Mouse Monoclonal anti-Mouse Qa-2 Antibody
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	C3H.SW mouse skin graft and splenocytes
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.