

Product Datasheet

Rat Monoclonal anti-Mouse Ly-6C Antibody (orb1153934)

Description	Rat Monoclonal anti-Mouse Ly-6C Antibody
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC
Immunogen	L3 cloned CTL cells
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rat IgG2c
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.