

Product Datasheet

Rat Monoclonal anti-Mouse Ly-6A/E (Sca-1) Antibody (orb1153937)

Description	Rat Monoclonal anti-Mouse Ly-6A/E (Sca-1) Antibody
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP, WB
Immunogen	IL-2-dependent mouse T-cell line (CTL-L)
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rat IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.