

Product Datasheet

Rat Monoclonal anti-Mouse IL-17A Antibody (orb1153964)

Description	Rat Monoclonal anti-Mouse IL-17A Antibody
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISpot, FC
Immunogen	E. coli expressed, recombinant mouse IL-17A
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rat IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.