

Product Datasheet

Mouse Recombinant anti-Mouse CD45 Recombinant Antibody (orb1154066)

Description	Mouse Recombinant anti-Mouse CD45 Recombinant Antibody
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP, WB
Immunogen	Mouse thymus or spleen
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze- thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Recombinant
Expiration Date	12 months from date of receipt.