

Product Datasheet

Rat Monoclonal anti-Mouse CD45 Antibody (orb1154067)

Description	Rat Monoclonal anti-Mouse CD45 Antibody
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP, WB
Immunogen	Mouse thymus or spleen
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rat IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.