

Product Datasheet

Rat Monoclonal anti-Human CD80 Antibody (orb1154376)

Description	Rat Monoclonal anti-Human CD80 Antibody
Species/Host	Rat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	RBL1 cells stably transfected with human CD80 (full length)
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rat IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.