

Product Datasheet

Mouse Recombinant anti-Human CD163 Recombinant Antibody (orb1154586)

Description	Mouse Recombinant anti-Human CD163 Recombinant Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC
Immunogen	Recombinant human CD163
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze- thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Recombinant
Expiration Date	12 months from date of receipt.