


Product Datasheet

VCAM1 antibody (orb11549)

Description	VCAM1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Mouse VCAM 1 protein
Target	VCAM1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	81 kDa
Uniprot ID	P19320
Dilution Range	WB=1:500-1000
Expiration Date	12 months from date of receipt.



IHC-P of human glomerulonephritis tissue...