



Product Datasheet

VEGF antibody (orb11553)


Description	VEGF antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human VEGF. Please contact us for the exact immunogen sequence. The peptide is available as orb374754.
Target	VEGF
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	24 kDa
Uniprot ID	P15692
NCBI	001020537, 21, 001025366
Entrez	7422
Dilution Range	IHC-P: 1:100-500 (based on 0.5 mg/ml), WB: 1:100-300
Expiration Date	12 months from date of receipt.



Western blot analysis of rat small intes...



IHC-P staining of mouse brain tissue usi...



IHC-P image of mouse liver tissue using ...