

Product Datasheet

VEGF antibody (orb11554)

Description

VEGF antibody

Species/Host

Rabbit

Conjugation

Unconjugated

Tested Applications

ELISA, IHC-P, WB

Immunogen

KLH conjugated synthetic peptide derived from human VEGF. Please contact us for the exact immunogen sequence. The peptide is available as orb374755.

Target

VEGF

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

23 kDa

Uniprot ID

[P15692](#), [Q00731](#), [P16612](#)

NCBI

[001020537](#), [21](#), [001025366](#)

Entrez



[2321](#)

Dilution Range


IHC-P: 1:100-800 (based on 0.5 mg/ml), WB: 1:100-500

Expiration Date


12 months from date of receipt.



IHC-P image of mouse liver tissue using ...



WB analysis of rat small intestines (lan...



Immunohistochemical staining of paraffin...