



Product Datasheet

VIP antibody (orb11560)


Description	VIP antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human VIP (125-152/170 aa)
Target	VIP
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	3.3/18 kDa
Uniprot ID	P01282
Dilution Range	IHC-P=1:100-500, IHC-F=1:100-500
Expiration Date	12 months from date of receipt.



orb11560 shows VIP expression in rat bre...



IHC-P of rat breast tissue (orb11560 at ...



Immunohistochemical staining of paraffin...