

Product Datasheet

Na⁺/K⁺-ATPase α 1 (phospho Ser16) Polyclonal Antibody (orb1156368)

Description	Na ⁺ /K ⁺ -ATPase α 1 (phospho Ser16) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.