

## Product Datasheet

### eEF2K (phospho Ser359) Polyclonal Antibody (orb1156875)

|                            |                                               |
|----------------------------|-----------------------------------------------|
| <b>Description</b>         | eEF2K (phospho Ser359)<br>Polyclonal Antibody |
| <b>Species/Host</b>        | Rabbit                                        |
| <b>Reactivity</b>          | Human, Mouse, Rat                             |
| <b>Conjugation</b>         | Unconjugated                                  |
| <b>Tested Applications</b> | ELISA, IHC-P                                  |
| <b>Note</b>                | For research use only                         |
| <b>Clonality</b>           | Polyclonal                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.               |