

## Product Datasheet

# HDAC6 (phospho Ser22) Polyclonal Antibody (orb1158949)

|                            |                                              |
|----------------------------|----------------------------------------------|
| <b>Description</b>         | HDAC6 (phospho Ser22)<br>Polyclonal Antibody |
| <b>Species/Host</b>        | Rabbit                                       |
| <b>Reactivity</b>          | Human, Mouse                                 |
| <b>Conjugation</b>         | Unconjugated                                 |
| <b>Tested Applications</b> | ELISA, IF, IHC-P, WB                         |
| <b>Note</b>                | For research use only                        |
| <b>Clonality</b>           | Polyclonal                                   |
| <b>Expiration Date</b>     | 12 months from date of receipt.              |