

Product Datasheet

GAPDH antibody (orb11605)

Description

Goat polyclonal to GAPDH (glyceraldehyde 3-

important en...

Species/Host

Goat

Reactivity

Bovine, Canine, Donkey, Equine, Feline, Fish, Gallus, Goat, Guinea pig, Hamster, Human, Mouse, Other, Porcine, Primate, Rabbit, Rat, Sheep

Conjugation

Unconjugated

Tested

IF, IHC-Fr, IHC-P, WB

Applications
Immunogen

Purified recombinant peptide derived from within residues 240 aa to the C-terminus of human GAPDH produced in E. coli.

Target

GAPDH

Preservatives

PBS, 20% glycerol and 0.05% sodium azide

Form/Appearance

Polyclonal antibody supplied as a 200 or 500 µl (2 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.

Concentration

2 mg/ml

Storage

Store at -20°C for long-term storage. Store at 2-8°C for up to one month.

Note

For research use only

Application notes

The antibody solution should be gently mixed before use.

Isotype

IgG

Clonality


Polyclonal

Dilution Range


WB:1:500-1:5,000, IF:1:50-1:250, IHC-P:1:200-1:1,000, IHC-F:1:200-1:1,000

Expiration Date

12 months from date of receipt.



Western blot analysis of 293 HEK cell ly...



Immunofluorescence analysis of Hepa1-6 C...