

## Product Datasheet

# LIMK-1/2 (phospho Thr508/505) Polyclonal Antibody (orb1160651)

|                            |                                                   |
|----------------------------|---------------------------------------------------|
| <b>Description</b>         | LIMK-1/2 (phospho Thr508/505) Polyclonal Antibody |
| <b>Species/Host</b>        | Rabbit                                            |
| <b>Reactivity</b>          | Human, Mouse, Rat                                 |
| <b>Conjugation</b>         | Unconjugated                                      |
| <b>Tested Applications</b> | ELISA, IHC-P, WB                                  |
| <b>Note</b>                | For research use only                             |
| <b>Clonality</b>           | Polyclonal                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                   |