


## Product Datasheet

### GST antibody (orb11611)

|                            |                                                                                                                                                                                |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Goat polyclonal antibody to GST (Glutathione-S-Transferase). GST is a 26 kDa protein encoded by the ...                                                                        |
| <b>Species/Host</b>        | Goat                                                                                                                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                   |
| <b>Tested Applications</b> | WB                                                                                                                                                                             |
| <b>Immunogen</b>           | Purified recombinant protein produced in E. coli.                                                                                                                              |
| <b>Target</b>              | GST                                                                                                                                                                            |
| <b>Preservatives</b>       | PBS, 20% glycerol and 0.05% sodium azide                                                                                                                                       |
| <b>Form/Appearance</b>     | Polyclonal antibody supplied as a 200 or 500 µl (6 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum. |
| <b>Concentration</b>       | 6 mg/ml                                                                                                                                                                        |
| <b>Storage</b>             | Store at -20°C for long-term storage. Store at 2-8°C for up to one month.                                                                                                      |
| <b>Note</b>                | For research use only                                                                                                                                                          |
| <b>Application notes</b>   | The antibody solution should be gently mixed before use.                                                                                                                       |
| <b>Isotype</b>             | IgG                                                                                                                                                                            |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                     |
| <b>Dilution Range</b>      | WB:1:500-1:2,000                                                                                                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                |



Western blot analysis of transformed cel...