



Product Datasheet

GFP antibody (orb11615)

Description	Goat polyclonal antibody to GFP (green fluorescent protein). GFP is a protein composed of 238 amino ...
Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, IHC-P, WB
Immunogen	Purified recombinant peptide produced in E. coli.
Target	In 293HEK cells transfected with cds plasmid detects a band of 27 kDa by Western blot. This antibody does not recognize mCherry fluorescent protein
Preservatives	PBS, 20% glycerol and 0.05% sodium azide
Form/Appearance	PBS, 20% glycerol and 0.05% sodium azide
Concentration	3 mg/ml
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only
Application notes	The antibody solution should be gently mixed before use.
Isotype	IgG
Clonality	Polyclonal
Purity	This antibody is epitope-affinity purified from goat antiserum
Dilution Range	WB:1:500-1:5,000, IF:1:50-1:1,000, IHC-P:1:50-1:1,000, IHC-F:1:50-1:1,000, IEM:1:50-1:1,000
Expiration Date	12 months from date of receipt.



Western blot analysis of transfected 293...



Confocal immunofluorescence analysis of ...