

Product Datasheet

Rab5a antibody (orb11628)

Description

Goat polyclonal antibody to mouse Rab5a. Rab5a

Species/Host

Goat

Reactivity

Bovine, Canine, Donkey, Equine, Feline, Gallus, Goat, Guinea pig, Hamster, Human, Mouse, Other, Porcine, Primate, Rabbit, Rat, Sheep

Conjugation

Unconjugated

Tested

IF, IHC-Fr, IHC-P, WB

Applications
Immunogen

Purified recombinant peptide derived from within residues 115 aa to the C-terminus of mouse Rab5a produced in E. coli.

Target

Rab5a

Preservatives

PBS, 20% glycerol and 0.05% sodium azide

Form/Appearance

Polyclonal antibody supplied as a 200 µl (4 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.

Concentration

4 mg/ml

Storage

Store at -20°C for long-term storage. Store at 2-8°C for up to one month.

Note

For research use only

Application notes

The antibody solution should be gently mixed before use.

Isotype

IgG

Clonality


Polyclonal

Dilution Range


WB:1:250-1:5,000, IF:1:50-1:250, IHC-P:1:100-1:500, IHC-F:1:100-1:500

Expiration Date

12 months from date of receipt.



Western blot analysis of transfected 293...



Confocal immunofluorescence analysis of ...