

## Product Datasheet

# Histone H3 (Tri Methyl Lys79) Monoclonal Antibody(3G3) (orb1162972)

|                            |                                                           |
|----------------------------|-----------------------------------------------------------|
| <b>Description</b>         | Histone H3 (Tri Methyl Lys79)<br>Monoclonal Antibody(3G3) |
| <b>Species/Host</b>        | Mouse                                                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                         |
| <b>Conjugation</b>         | Unconjugated                                              |
| <b>Tested Applications</b> | IF, IHC-P, IP, WB                                         |
| <b>Note</b>                | For research use only                                     |
| <b>Clonality</b>           | Monoclonal                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                           |