

Product Datasheet

Rab25 antibody (orb11633)

Description

RAB25 belongs to the large RAB family of low

Species/Host

Goat

Reactivity

Bovine, Canine, Donkey, Equine, Feline, Gallus, Goat, Guinea pig, Hamster, Human, Mouse, Other, Porcine, Primate, Rabbit, Rat, Sheep

Conjugation

Unconjugated

Tested

IF, IHC-Fr, IHC-P, WB

Applications
Immunogen

Purified recombinant peptide derived from within residues 120 aa to the C-terminus of mouse Rab25 produced in E. coli.

Target

Rab25

Preservatives

PBS, 20% glycerol and 0.05% sodium azide

Form/Appearance

Polyclonal antibody supplied as a 200 µl (2 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.

Concentration

2 mg/ml

Storage

Store at -20°C for long-term storage. Store at 2-8°C for up to one month.

Note

For research use only

Application notes

The antibody solution should be gently mixed before use.

Isotype

IgG

Clonality


Polyclonal

Dilution Range


WB:1:250-1:2,000, IF:1:25-1:250, IHC-P:1:200-1:1,000, IHC-F:1:200-1:1,000

Expiration Date

12 months from date of receipt.



Western blot analysis of transfected 293...



Immunofluorescence analysis of Hepa1-6 C...