

## Product Datasheet

### N/H/K-Ras Polyclonal Antibody (orb1163354)

|                            |                                 |
|----------------------------|---------------------------------|
| <b>Description</b>         | N/H/K-Ras Polyclonal Antibody   |
| <b>Species/Host</b>        | Rabbit                          |
| <b>Reactivity</b>          | Human, Mouse, Rat               |
| <b>Conjugation</b>         | Unconjugated                    |
| <b>Tested Applications</b> | ELISA, IHC, WB                  |
| <b>Note</b>                | For research use only           |
| <b>Clonality</b>           | Polyclonal                      |
| <b>Expiration Date</b>     | 12 months from date of receipt. |