



Product Datasheet

AFP antibody (orb11644)


Description	Rabbit polyclonal antibody to AFP
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human Alpha-Fetoprotein (AFP)(NP_001125.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	68kDa
Uniprot ID	P02771
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cells...



Immunofluorescence analysis of human liver...



Western blot analysis of extracts of mouse liver...