

## Product Datasheet

### DNAM-1 (Phospho-Ser329) Polyclonal Antibody (orb1164506)

|                            |                                                |
|----------------------------|------------------------------------------------|
| <b>Description</b>         | DNAM-1 (Phospho-Ser329)<br>Polyclonal Antibody |
| <b>Species/Host</b>        | Rabbit                                         |
| <b>Reactivity</b>          | Human, Monkey, Mouse                           |
| <b>Conjugation</b>         | Unconjugated                                   |
| <b>Tested Applications</b> | ELISA, IF, IHC, WB                             |
| <b>Note</b>                | For research use only                          |
| <b>Clonality</b>           | Polyclonal                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                |