

## Product Datasheet

### FAK (Phospho-Tyr576) Antibody (orb1164868)

|                            |                                  |
|----------------------------|----------------------------------|
| <b>Description</b>         | FAK (Phospho-Tyr576)<br>Antibody |
| <b>Species/Host</b>        | Rabbit                           |
| <b>Reactivity</b>          | Human, Mouse, Rat                |
| <b>Conjugation</b>         | Unconjugated                     |
| <b>Tested Applications</b> | WB                               |
| <b>Note</b>                | For research use only            |
| <b>Clonality</b>           | Polyclonal                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |